




# ALLERG-O-SYSTEM

## Manual test procedure



### Features for *in-vitro* allergy diagnostics

- ALLERG-O-SYSTEM including reader, washer, light table (optional)
- Quantitative determination of specific IgE, IgG, IgG4 and Total IgE
- 4 or 6 point calibration (WHO 75/502 calibrated) for specific IgE
- Up to 800 determinations per day
- Optimized for the use with panels

### Related Products for allergy *in-vitro* diagnostics

- |  |                        |                          |                      |
|--|------------------------|--------------------------|----------------------|
| • Specific IgE EAST with 4 Calibrators | <b>REF</b> 0540200PKL& | • Specific IgE REAST     | <b>REF</b> 052xxxxFL |
| • Reference System                     | <b>REF</b> 0541000PKL  | • Total IgE ALLERG-O-LIQ | <b>REF</b> 08101FL   |
| • Specific IgE EAST with 6 Calibrators | <b>REF</b> 074000PQ    | • Reference System       | <b>REF</b> 07050FL   |
| • Reference System                     | <b>REF</b> 0560200PKL& | • Allergens              | <b>REF</b> Code-FL   |
| • Allergens                            | <b>REF</b> 0561000PKL  | • Specific IgG           | <b>REF</b> 10100PG   |
| • Total IgE-HRP EIA                    | <b>REF</b> 076000PQ    | • Specific IgG4          | <b>REF</b> 10104PG   |
|  | <b>REF</b> Code        | • Allergens              | <b>REF</b> 13-Code-G |
|  | <b>REF</b> 08102CP     |                          |                      |

**DR. FOOKE**  
LABORATORIEN GmbH

2011-05